WHOMS2 Software Installation Guide

Step 1

• Copy cable driver and WHOMS2 to C drive app folder.

Step 2

• Install the USB To Serial Driver.

Step 3

• Locate WHOMS2 folder and click on WHOMS2 Setup.exe

Name	Date modified	Туре	Size
D tech driver	8/5/2024 10:10 AM	File folder	
PL2303_Prolific_DriverInstaller_v1.12.0	8/12/2024 11:43 AM	File folder	
WHOMS2	12/8/2023 1:35 PM	File folder	
🚮 0x0409	9/12/2023 11:08 AM	Configuration setti	22 KB
how to install WHOMS2 in laptop	8/5/2024 3:01 PM	Microsoft Word D	165 KB
🐸 setup	9/12/2023 11:08 AM	Application	1,263 KB
📲 Setup	9/12/2023 11:08 AM	Configuration setti	6 KB
whoms D Tech cable driver	8/20/2024 1:00 PM	Microsoft Word D	1,511 KB
🛃 WHOMS2 Setup	9/12/2023 11:08 AM	Windows Installer	1,422 KB

Step 4

• Click the **Next** button.



Step 5

• Click the Change button



Step 6

Browse the file location (c:\Program Files (x86)\WHOMS2\)



Step 7

Click Next

	ale al an an a constant							
🛃 WHOMS	2 Setup - InstallShield Wizar	d		×				
Destination Click Nex	on Folder at to install to this folder, or click t	Change to install to	hange to install to a different folder.					
	Install WHOMS2 Setup to: C:\Program Files (x86)\WHOM:	S2\		Change				
InstallShield -								
		< Back	Next >	Cancel				

Step 8

Click Install and you're done installing.

JWHOMS2 Setup - InstallShield Wizard	X
Ready to Install the Program The wizard is ready to begin installation.	5
If you want to review or change any of your installation settings, click Back. Click Cancel to exit the wizard. Current Settings:	
Setup Type: Typical	
Destination Folder: C:\Program Files (x86)\WHOMS2\	
User Information: Name: user Company:	
InstallShield Cancel	

WHOMS2 Software Guide

Before launching WHOMS2 software, kindly make sure that you have these following things.

CSV Folder

The CSV folder is the folder that contains all the data and measurements of the students. It is therefore very important to handle it with precautions. Please take some time to read the following advice regarding the CSV folder.

1. At all times, make sure the folder contains **ONLY the Height & Weight** CSV Files.

2. Avoid opening the folder directly (via Explorer) to minimize corruption of any file.

3. If need to backup any file or the CSV Folder, do not save it within the CSV folder.

Connecting B1000-S BMI Machine

If the machine is properly **connected**, you will see a **green light** on the top right side of the WHOMS window.



If you see only **red light**, please close the WHOMS program, **check the connection** and restart WHOMS. If all attempts fail, please confirm the driver of the cable has been properly installed.

🐳 WHOMS2 - Canto	onment Primary School					- 6	×
Open F	older Generate Report	Generate TAF	AnBbCol AnBo of A 1110000 (1100 Help	AaBbt AaBbC Ieiding 1 Heading 2 Ieiding Exit	AaBbCc AaBbCc Heading 3 Heading	Replace	•
No.	Name	Gender Age	Status Height	Weight			
					Height		
					Weight		
					A	Commit	
					<	Skip	
						Discard	

Guide

NHOMS2 - Cantonment Prin	mary School						- 0	×
Open Folder	Generate Report	Generate TAF	He	lp	Exit			
No.	Name	Gender Ar	ne Status X	Height	Weight			
		Select the folder that contains the C	SV files			Hei	ght	
		Desktop A Home Gallery OneDrive - Personal Desktop Desktop Documents Cogniles Chate created: 9/9/20 PDFX	124 2:14 PM			Wei	ight	
						Α	Commit	
						<	Skip	
					1		Discard	

1. Click "**Open Folder**" and select the folder with the extracted CSV Files.

2. Select a file from the drop-down list (grey bar) to conduct measurement.

WHOMS2 - Cantonment Primary School		- 🗆 X
Open Folder Generate Report Generate TAF Help Exit	1	•
P1-M_Height_Weight_20240628070246.csv 🔉		
P6-C_Height_Weight_20240628070304.csv	1	
	Hei	ght
	Wei	ight
	A	Commit
	<	Skip
		Discard

3. The class list should be loaded into the data area. If you don't see a "**1**" at the registration no. box, key in manually.

WHC	JMS2 - Cantonment Primary School								x		
	Open Folder Generate Report (Generate T/	MF	He	elp	Б	xit		Registration no. of the student		
P1-M	_Height_Weight_20240628070246.csv										
No.	Name	Gender	Age	Status	Height	Weight		ANG JIA YING, GRACE			
1	ANG JIA YING, GRACE	F	6	м	165	64.3					
2	AURINA MA HUI XIN	F	7	м				Hei	ght		
3	CHLOE LIM XINYU	F	7	м					ight		
4	CHUA XIN TONG	F	7	м				vve	igni		
5	FANG MINXIN GLORIA	F	6	м							
6	FANG MINZHI NORA	F	6	м							
7	FONG JENG MUN VICTORIA	F	7	м							
8	HAZEL LEE EN SHU	F	6	м							
9	ONG NING ZHEN GABRIELLE	F	6	м							
10	QUEK YING LENG, AUDREY	F	6	м							
11	STEPHY SOPHIA NG CHU YUN	F	6	м							
12	TAN XIAN QING	F	7	м				A	Commit		
13	TAN YI ROU FAITH	F	6	м				<	Skip		
14	TERESA EE JING QI	F	6	м					Sub		
15	ANG LI JIE LUKE	М	7	м					Discard		

5. Ask the first student to stand on the machine. Ask the student to step down when the incoming data is acceptable.

WHO	0MS2 - Cantonment Primary School								- 🗆 X	
P1-M	Open Folder Generate Report	Generate T/	ΛF	He	əlp	Exi	t		2	
No.	No. Name Gender Age Status Height Weight							AURINA MA HUI XIN		
1	ANG JIA YING, GRACE	F	6	м	165	64.3				
2	AURINA MA HUI XIN	F	7	M				Height	165	
3	CHLOE LIM XINYU	F	7	м				Maight	(0.7	
4	CHUA XIN TONG	F	7	м				weight	63./	
5	FANG MINXIN GLORIA	F	6	м				Severely Overweight		
6	FANG MINZHI NORA	F	6	м						
7	FONG JENG MUN VICTORIA	F	7	м				BMI: 23.4		
8	HAZEL LEE EN SHU	F	6	м			1	Ideal: 35.9	kg - 52kg	
9	ONG NING ZHEN GABRIELLE	F	6	м				-11.7kg to	Acceptable	
10	QUEK YING LENG, AUDREY	F	6	м				4 1kg to O	vorwoight	
11	STEPHY SOPHIA NG CHU YUN	F	6	м				-4. TKY 10 O	verweight	
12	TAN XIAN QING	F	7	м				A	Commit	
13	TAN YI ROU FAITH	F	6	м				<	Skip	
14	TERESA EE JING QI	F	6	м					-mb	
15	ANG LI JIE LUKE	М	7	м					Discard	

7. Click "Commit" to commit the measured data.

WHC	🐳 WHOM52 - Cantonment Primary School — 🗆 🗙										
	Open Folder Generate Report	Generate T <i>A</i>	ΛF	He	эlр	Exit		2			
P1-M	_Height_Weight_20240628070246.csv	·									
No.	Name	Gender	Age	Status	Height	Weight	AURINA	A MA HUI XIN			
1	ANG JIA YING, GRACE	F	6	М	165	64.3					
2	AURINA MA HUI XIN	F	7	м			Height	165			
3	CHLOE LIM XINYU	F	7	М			Woight	10.7			
4	CHUA XIN TONG	F	7	М			weight	63./			
5	FANG MINXIN GLORIA	F	6	М			Severely	Overweight			
6	FANG MINZHI NORA	F	6	М			DMI: 02.4	overweight			
7	FONG JENG MUN VICTORIA	F	7	М			BIVII: 23.4				
8	HAZEL LEE EN SHU	F	6	М			Ideal: 35.9	kg - 52kg			
9	ONG NING ZHEN GABRIELLE	F	6	М			-11.7kg to	Acceptable			
10	QUEK YING LENG, AUDREY	F	6	М			-1 1kg to C)verweight			
11	STEPHY SOPHIA NG CHU YUN	F	6	М			-4. IKY to C	verweight			
12	TAN XIAN QING	F	7	М			A	Ççmmit			
13	TAN YI ROU FAITH	F	6	М				Commit the data to the resp			
14	TERESA EE JING QI	F	6	М							
15	ANG LI JIE LUKE	М	7	М				Discard			

You can also measure a specific student

- 1. Enter the registration number directly into the Registration Number Box.
- 2. Click "Skip" to move to the next student and "<" to the previous student.
- 3. Click directly on the student in the class list.

Save data

WHOMS automatically saves the measured data every time the class is being **CHANGED**. It also saves the data when the program is closed.

Hence, please do not allow the program to terminate abruptly; this may result in loss of data.

Report Generation

Report Generation works on the currently opened folder. If the current folder only has a level of student's data, report generated would only be based on that level

Once you have done with measurements, you can click on "Generate Report" option.

	DMS2 - Cantonment Primary School								- 🗆	×
P1-M	Open Folder Generate Report Height_Weight_20240	Senerate TA	ιF	He	эр	Exit	~		4	•
No.	Name	Gender	Age	Status	Height	Weight	1	CHUA	XIN TONG	
1	ANG JIA YING, GRACE	F	6	м	165	64.3				
2	AURINA MA HUI XIN	F	7	м	165	63.7	H	leight		
3	CHLOE LIM XINYU	F	7	м	165	63.8	1 10	laight		
4	CHUA XIN TONG	F	7	м				veigni		
5	FANG MINXIN GLORIA	F	6	м						
6	FANG MINZHI NORA	F	6	м						
7	FONG JENG MUN VICTORIA	F	7	м						
8	HAZEL LEE EN SHU	F	6	м			'			
9	ONG NING ZHEN GABRIELLE	F	6	м						
10	QUEK YING LENG, AUDREY	F	6	м						
11	STEPHY SOPHIA NG CHU YUN	F	6	м						
12	TAN XIAN QING	F	7	м			A		Commit	
13	TAN YI ROU FAITH	F	6	м			<		Skip	
14	TERESA EE JING QI	F	6	м					Sinp	
15	ANG LI JIE LUKE	М	7	М					Discard	

Your report will generate and will automatically save onto the user's desktop folder with the naming convention "Report ddMMyyyy HHmmss".



🗴 AutoSave 🔵 off) 🖫 🍤 - 🖓 - 🖉 Report 10092024110221 • Saved to this PC												Search
I	-ile	Home	Insert	Page Layout	Formulas	Data	Review	Vie	w Au	ıtomate	Help	
ſ	ĥ	X	Arial	~8	~ A^ A*	$\equiv \Xi$	≡ = ²	% ~	ab c	Wrap Tex	t	General
	Paste	L <u>₽</u> ~ ≪	B I	<u>u</u> • = •	<u>∽</u> ~ <u>A</u> ~			-= -=		Merge &	Center 🗸	\$ ~ 9
	Clipbo	ard 🗔		Font	۲	i		Ali <u>c</u>	gnment		I	151 N
	D SEC	CURITY W	ARNING	Some active conte	ent has been dis	abled. C	lick for mo	ore det	ails.	Enable Co	ntent	
J1		~	: X ~	fx								
		Α	B		С		D	E	F	G	н	J
1		Class	Reg No.		Name		Gender	Age	Status	Height	Weight	
2	P1-M		1	ANG JIA YING, GI	RACE		F	6	М	165	64.3	÷
3	P1-M		2	AURINA MA HUI X	an		F	7	М	165	63.7	
4	P1-M		3	CHLOE LIM XINYU	J		F	7	М	165	63.8	

Open the 'XLSM' file and navigate to 'Raw Data' tab to display saved data.

More questions

Further enquiries, please email to itd@avamech.com or call 6746 8062 during Monday to Friday, 0900 – 1700.